

PAR PULSE LAVAGE SYSTEM



- With new design PAR pulse lavage system provides a wider viewing, provides effective and quick cleaning of bone pieces necrotic tissues, bacteria and foreign items.
- It has a longer lasting performance.
- Pulse Lavage It is completely disposable.




- Ergonomic design comfortable to handle and use.
- Battery pack easy to remove user friendly design.
- The two-stage speed setting meets all special needs.
- Battery built-in.




-The pulse lavage has the ability to significantly clean and remove debris, whilst simultaneously applying suction.

-Easy attached & released tip accessor.

- Different tips for knee and hip are available.



- Dual spike & Single options are available.



Ortiz Medikal Tur. İnş. San. Ve Tic. Ltd. Şti.
Doğanlar Mah. 1411/1 Sok. No:28/A
Bornova/İzmir/Turkey

www.ortizmedical.com.tr

T: +90 232 435 5546

F: +90 232 435 5517

hale@ortizmedical.com.tr